

# Southern Group Laboratory

## SAFETY DATA SHEET

### 1. Identification of the substance/preparation and of the company/undertaking

#### *Identification of the product*

SGL Product Code: 3244  
Product Name: Fluorocult Lauryl Sulphate broth

#### *Manufacturer/supplier identification*

Company Name: SGL Ltd., Cavendish Courtyard, Sallow Road, Weldon North Ind. Estate,  
Corby, Northants. NN17 5JX  
Email: info@sglab.co.uk  
Emergency Telephone No.: (01536) 403815 Available Monday to Friday 9-5

### 2. Composition/Information on Ingredients

This medium does not contain any substances presenting a health hazard within the meaning of the Dangerous Substances Directive

### 3. Hazards Identification

Not classified as dangerous according to EC Directives

### 4. First Aid Measures

Eye contact: Irrigate thoroughly with plenty of water.  
Skin contact: Wash off skin thoroughly with soap and water.  
Ingestion: Wash out mouth thoroughly with water. In severe cases, obtain medical attention.

### 5. Firefighting Measures

#### *Special risks:*

Not combustible

### 6. Accidental Release Measures

Wear appropriate protective clothing.  
Transfer into suitable containers for recovery or disposal. Flush area with plenty of water.

### 7. Handling and Storage

#### *Handling:*

Avoid contact with eyes, skin and clothing

#### *Storage:*

Store at 2 - 25°C

## **8. Exposure Controls / Personal Protection**

### ***UK exposure limits:***

None assigned

### ***Personal protective equipment***

Gloves: Disposable vinyl

## **9. Physical and Chemical Properties**

### ***General information:***

Form: Liquid  
Colour: Yellowish/Brown  
Odour: Odourless

## **10. Stability and Reactivity**

Stable under normal conditions

Avoid exposure to direct sunlight

## **11. Toxicological Information**

Low order of acute toxicity

## **12. Ecological Information**

No relevant studies identified.

## **13. Disposal Considerations**

Sterilise used containers by autoclaving before final disposal. Dispose of in accordance with applicable local and national regulations.

## **14. Transport Information**

Not subject to transport regulations

## **15. Regulatory Information**

### ***Labelling according to EC directives***

Not classified as dangerous according to EC Directives.

## **16. Other Information**

Issue Date: 26/10/2015  
Issue No: 01