

SOUTHERN GROUP LABORATORY LIMITED

SAFETY DATA SHEET

SECTION 1.	IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING
-------------------	---

1.1 Product identifier

Product name Chromogenic Listeria Selective Agar

SGL Product code 8126

This product is a mixture

7447.41 Relevant identified uses of the substance or mixture and uses advised against

Identified uses Reagent for laboratory analysis

Uses advised against No information on uses against available

1.3. Details of the supplier of the safety data sheet

Supplier Southern Group Laboratory,
Cavendish Courtyard, Sallow Road,
Weldon North Industrial Estate,
Corby, Northants,
United Kingdom
NN17 5JX

TEL: +44(0) 1536 403815

FAX: +44(0) 1536 403814

EMAIL: info@sglab.co.uk

Available 9am-5pm (UK time) Mon to Fri except UK public holidays

Emergency telephone NHS Direct: 0845 4647 (England and Wales)
NHS 24: 08454 24 24 24 (Scotland)

SECTION 3.	COMPOSITION/INFORMATION ON INGREDIENTS
-------------------	---

3.1. Mixtures

Name	EC No.	CAS-No.	Concentration	Classification
Lithium chloride	231-212-3	7447-41-8	<2% w/v	Acute Toxicity 4 (oral):H302 Skin Irritation 2: H315 Eye Irritation 2: H319
Cycloheximide	200-636-0	66-81-9	<0.01% w/v	Germ Cell Mutagenicity 2: H341 Reproductive Toxicity 1B: H360D Acute Toxicity 2 (oral):H300 Hazardous to the Aquatic Environment, Chronic 2: H411
Nalidixic acid	206-864-7	389-08-2	<0.01% w/v	Respiratory Sensitisation 1: H334 Skin Sensitisation 1: H317
Polymyxin B sulphate	215-774-7	1405-20-5	<0.01% w/v	Acute Toxicity 4 (oral):H302
Ceftazidime hydrate	616-626-8	78439-06-2	<0.01% w/v	Respiratory Sensitisation 1: H334 Skin Sensitisation 1: H317
Sodium carbonate	207-838-8	497-19-8	<4% w/v	Eye Irritation 2: H319

Composition comments The data shown are in accordance with Regulation (EC) No. 1272/2008.

SECTION 4.	FIRST-AID MEASURES
-------------------	---------------------------

4.1 Description of first aid measures

Inhalation Remove from exposure

Ingestion Wash out mouth thoroughly with water. In severe cases and/or if symptoms occur, obtain medical attention.

Skin contact Wash off thoroughly with soap and water. Remove contaminated clothing and wash before re-use.

Eye Contact Irrigate thoroughly with water for at least 15 minutes. Obtain medical attention if symptoms persist.

Self Protection of the first aider No special precautions identified

4.2 Most important symptoms and effects, both acute and delayed

None expected under normal circumstances. In extreme cases may cause nausea, vomiting and/or diarrhoea if ingested.

4.3 Indication of any immediate attention and special treatment needed

Notes for doctor No specific recommendations

SECTION 5.	FIRE – FIGHTING MEASURES
-------------------	---------------------------------

5.1. Extinguishing media

Suitable extinguishing media Water fog, foam, carbon dioxide (CO₂), dry powder

Unsuitable extinguishing media None identified

Prevent fire fighting media from entering watercourses or ground water.

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products Possible generation of hazardous gases or vapours in event of fire

5.3. Advice for firefighters

Protective actions during firefighting Wear self contained breathing apparatus (SCBA)

SECTION 6.	ACCIDENTAL RELEASE MEASURES
-------------------	------------------------------------

6.1. Personal precautions, protective equipment and emergency procedures

Wear personal protective clothing as described in Section 8 of this safety data sheet.

6.2. Environmental precautions

Should not be released into the environment

6.3. Methods and material for containment and cleaning up

Collect and seal material in closed container. Wear personal protective clothing as recommended in Section 8. Wash site of spillage thoroughly with water and detergent. Dispose as recommended in Section 13.

6.4. Reference to other sections

Section 8 for protective measures
Section 13 for disposal considerations

SECTION 7.	HANDLING AND STORAGE
-------------------	-----------------------------

7.1 Precautions for safe handling

Wear personal protective clothing. Avoid contact with skin, eyes and clothing. Change contaminated clothing. Wash hands after working with material.

7.2 Conditions for safe storage, including any incompatibilities

Store according to label requirements. Keep closed.

7.3 Specific end use(s)

Use in laboratories

SECTION 8.	EXPOSURE CONTROLS/PERSONAL PROTECTION
-------------------	--

8.1 Control parameters

Component	CAS No	Long Term Exposure Limit (8 hr TWA Reference Period)	Short Term Exposure (15 Min Reference Period)	Basis	Comments
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable	Not applicable

This mixture contains no components with known occupational exposure limits.

8.2 Exposure controls

Engineering measures Avoid contact with skin, eyes and clothing. Wash hands before breaks.

Protective equipment



Eye/face protection None needed during normal use

Hand protection Disposable nitrile rubber gloves complying with EC Directive 89/686/EEC and Standard EN374
 Glove thickness of 0.11mm
 See manufacturer's recommendations for break through times
 Material tested: Not applicable

Other skin and body protection Long sleeved clothing

Hygiene measures Wash hands after use

Respiratory protection None needed during normal use

SECTION 9.	PHYSICAL AND CHEMICAL PROPERTIES
-------------------	---

9.1. Information on basic physical and chemical properties

Appearance	Solid
Colour	Cream / Yellow
Odour	No data available
Odour threshold	No data available
pH	7.0 to 7.4 at 20°C
Melting point/freezing point	No data available
Initial boiling point and boiling range	No data available
Flash point	No data available
Evaporation rate	No data available
Flammability (solid, gas)	No data available
Upper/lower flammability or explosive limits	
Lower explosion limit (VOL-%)	No data available
Upper explosion limit (VOL-%)	No data available
Vapour pressure	No data available
Vapour density	No data available
Relative density	No data available
Solubility(ies)	No data available
Partition coefficient: n-octanol/water	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Explosive properties	No data available
Oxidising properties	No data available

9.2. Other information

None

SECTION 10.	STABILITY AND REACTIVITY
--------------------	---------------------------------

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under storage at normal conditions

10.3 Possibility of hazardous reactions

None known

10.4 Conditions to avoid

None known

10.5 Incompatible materials

None known

10.6 Hazardous decomposition products

None known

SECTION 11. TOXICOLOGICAL INFORMATION
--

11.1 Information on toxicological effects

The acute toxicity data is based on the mixture.

Acute oral toxicity	Classification criteria are not met.
ATE_{mix}	>2000 mg/kg
Method	Calculation
Acute dermal toxicity	No data available
Acute inhalation toxicity	No data available
Skin corrosion/irritation	No data available
Serious eye damage/irritation	No data available
Respiratory or skin sensitisation	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
STOT single exposure	No data available
STOT repeated exposure	No data available
Aspiration hazard	No data available

11.2 Additional information

In extreme cases may cause nausea, vomiting and/or diarrhoea if ingested.

SECTION 12. ECOLOGICAL INFORMATION

12.1 Toxicity

No information available

12.2 Persistence and degradability

No information available

12.3 Bioaccumulative potential

No information available

12.4 Mobility in soil

No information available

12.5 Results of PBT and vPvB assessment

No information available

12.6 Other adverse effects

None known

SECTION 13.	DISPOSAL CONSIDERATIONS
--------------------	--------------------------------

13.1 Waste treatment methods

Sterilise used containers by, for example, autoclaving before final disposal. Dispose of in accordance with applicable local and national regulations.

Waste should not be disposed of by release to sewers.

SECTION 14.	TRANSPORT INFORMATION
--------------------	------------------------------

Not classified as dangerous according to transport regulations.

Land transport (ADR/RID) / Sea transport (IMDG) / Air transport (ICAO-TI/IATA-DGR)

14.1 UN number Not applicable

14.2 UN proper shipping name Not applicable

14.3 Transport hazard class(es) Not applicable

14.4 Packing group Not applicable

14.5 Environmental hazards No data available

14.6 Special precautions for user No data available

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable, packaged goods.

SECTION 15.	REGULATORY INFORMATION
--------------------	-------------------------------

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

EU legislation Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of Substances and mixtures (as amended).

15.2 Chemical safety assessment

No Chemical Safety Assessment has been carried out for this substance/mixture by the supplier.

SECTION 16.	OTHER INFORMATION
--------------------	--------------------------

16.1 Relevant H- and EUH-phrases (number and full text)

H300: Fatal if swallowed
H302: Harmful if swallowed
H315: Causes skin irritation
H317: May cause an allergic skin reaction
H319: Causes serious eye irritation
H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled
H360D: May damage fertility or the unborn child
H411: Toxic to aquatic life with long lasting effects

16.2 Additional information

The above information is based on present knowledge, but does not purport to be inclusive and should only be used as a guide. Southern Group Laboratory shall not be held liable for any damage resulting from handling or from contact with the above product.

The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made up material.

16.3 Revision summary

Changes to this issue Updated to format compliant with Regulation (EC) No 1907/2006 (REACH), the Globally Harmonized System of Classification and Labelling of Chemicals (GHS) and Regulation (EC) No 1272/2008 (CLP).

Revision Date: 02.03.2020
Issue No: 02

END OF SAFETY DATA SHEET