

SOUTHERN GROUP LABORATORY LTD

HEALTH AND SAFETY INFORMATION

Product:	Pseudomonas Selective CN Agar	Hazard Class	: NR
		UN No.	: n/a
		CAS No.	: n/a
SGL Code:	8383CN, 9095CN, 8096		

DESCRIPTION: Agar mixture with antibiotics

FIRE & EXPLOSION HAZARD: Not applicable

HEALTH HAZARD: May be harmful if ingested in quantity.
May irritate eyes.

Toxicity data: no data
Carcinogenicity: no evidence of effects
Mutagenicity/Teratogenicity: no evidence of effects
Exposure limits: not assigned

FIRST AID:

Eyes: Irrigate thoroughly with water for at least 10 minutes. If discomfort persists obtain medical attention.

Lungs: Remove from exposure, rest and keep warm. In severe cases obtain medical attention.

Skin: Wash off skin thoroughly with water. Remove contaminated clothing and wash before re-use. In severe cases obtain medical attention.

Mouth: Wash out mouth thoroughly with water, give plenty of water to drink. Obtain medical attention.

REACTIVE HAZARDS: none

ACCIDENTAL RELEASE: Wear appropriate protective clothing. Carefully sweep up and dispose of in accordance with local regulations. Wash site of spillage thoroughly with detergent and water.

PROTECTIVE MEASURES: (as appropriate to quantity handled)

Respirator: Not applicable
Ventilation: Not applicable
Gloves: Rubber or plastic
Eye Protection: Goggles or face shield
Other measures: Not applicable

STORAGE AND HANDLING: No special requirements