

# Southern Group Laboratory

## SAFETY DATA SHEET

### 1. Identification of the substance/preparation and of the company/undertaking

#### *Identification of the product*

SGL Product Code: 8390

Product name: Slanetz & Bartley Agar

#### *Manufacturer/supplier identification*

Company: SGL Ltd, Cavendish Courtyard, Sallow Road, Weldon North Ind. Est., Corby, Northants. NN17 5JX

Emergency telephone No: (01536) 403815

### 2. Composition/information on ingredients

This medium does not contain any substances presenting a health hazard within the meaning of the Dangerous Substances Directive

### 3. Hazards identification

Not classified as dangerous according to EC Directives

### 4. First aid measures

Eye contact: Irrigate thoroughly with plenty of water  
Skin contact: Wash off skin thoroughly with soap and water  
Ingestion: Wash out mouth thoroughly with water. In severe cases obtain medical attention.

### 5. Fire fighting measures

#### *Special risks:*

Not combustible

### 6. Accidental release measures

Wear appropriate protective clothing  
Transfer into suitable containers for recovery or disposal. Flush area with plenty of water

## **7. Handling and storage**

### ***Handling:***

Avoid contact with eyes, skin and clothing

### ***Storage:***

Store at 2 - 25°C

## **8. Exposure controls/personal protection**

### ***UK exposure limits:***

None assigned

### ***Personal protective equipment***

Gloves: Disposable vinyl

## **9. Physical and chemical properties**

### ***General information:***

Form:	Agar
Colour:	Clear, colourless or slightly pink
Odour:	Odourless

## **10. Stability and reactivity**

Stable under normal conditions

Avoid exposure to direct sunlight

## **11. Toxicological information**

Low order of acute toxicity

## **12. Ecological information**

No relevant studies identified

## **13. Disposal considerations**

Sterilise used containers by autoclaving before final disposal. Dispose of in accordance with applicable local and national regulations

## **14. Transport information**

Not subject to transport regulations

## **15. Regulatory information**

### *Labelling according to EC directives*

Not classified as dangerous according to EC Directives

## **16. Other information**

Issue Date: 11/10/2005

Date of print: 10/04/2008