

**SOUTHERN GROUP LABORATORY LTD
HEALTH AND SAFETY INFORMATION**

Product: Defibrinated Horse Blood	Hazard Class: NR
UN No.: n/a	
CAS No.: n/a	SGL Code: 9011

DESCRIPTION: 100% Horse Blood derived from healthy donor animals held in closed herds. Animals undergo regular veterinary inspection and are certified healthy.

FIRE & EXPLOSION HAZARD: Not applicable

HEALTH HAZARD: Non Hazardous. No significant health hazard is expected. To the best of our knowledge no hazard has been demonstrated with this material

Handle to minimise exposure for protection against unknown hazards

Possible routes of entry: Inhalation - Yes
Skin - Yes
Ingestion - Yes

Toxicity data: no data
Carcinogenicity: no evidence of effects
Mutagenicity/Teratogenicity: no evidence of effects
Exposure limits: not assigned

FIRST AID:

Eyes: Irrigate thoroughly with water for at least 10 minutes. If discomfort persists obtain medical attention.

Lungs: Remove from exposure, rest and keep warm. In severe cases obtain medical attention.

Skin: Wash off skin thoroughly with water.

Mouth: Wash out mouth thoroughly with water, give plenty of water to drink. OBTAIN MEDICAL ATTENTION.

REACTIVE HAZARDS: none

ACCIDENTAL RELEASE: Wear appropriate protective clothing. Carefully sweep up and dispose of in accordance with local regulations. Wash site of spillage thoroughly with detergent and water.

PROTECTIVE MEASURES: (as appropriate to quantity handled)

Respirator: Not applicable
Ventilation: Not applicable
Gloves: Rubber or plastic
Eye Protection: Goggles or face shield
Other measures: Not applicable

STORAGE AND HANDLING: No special measures required
Storage temperature.
On receipt store at 2 - 8°C

PHYSICAL AND CHEMICAL: Stable liquid containing horse plasma and red blood cells