

Southern Group Laboratory

SAFETY DATA SHEET

1. Identification of the substance/preparation and of the company/undertaking

Identification of the product

SGL Product Code: 9098

Product name: Slanetz & Bartley Agar

Manufacturer/supplier identification

Company: SGL Ltd, Cavendish Courtyard, Sallow Road, Weldon North Ind. Est., Corby,
Northants. NN17 5JX

Emergency telephone No: (01536) 403815

2. Composition/information on ingredients

This medium does not contain any substances presenting a health hazard within the meaning of the Dangerous Substances Directive

3. Hazards identification

Not classified as dangerous according to EC Directives

4. First aid measures

Eye contact: Irrigate thoroughly with plenty of water
Skin contact: Wash off skin thoroughly with soap and water
Ingestion: Wash out mouth thoroughly with water. In severe cases obtain medical attention.

5. Fire fighting measures

Special risks:

Not combustible

6. Accidental release measures

Wear appropriate protective clothing
Transfer into suitable containers for recovery or disposal. Flush area with plenty of water

7. Handling and storage

Handling:

Avoid contact with eyes, skin and clothing

Storage:

Store at 2 - 25°C

8. Exposure controls/personal protection

UK exposure limits:

None assigned

Personal protective equipment

Gloves: Disposable vinyl

9. Physical and chemical properties

General information:

Form:	Agar
Colour:	Clear, colourless or slightly pink
Odour:	Odourless

10. Stability and reactivity

Stable under normal conditions

Avoid exposure to direct sunlight

11. Toxicological information

Low order of acute toxicity

12. Ecological information

No relevant studies identified

13. Disposal considerations

Sterilise used containers by autoclaving before final disposal. Dispose of in accordance with applicable local and national regulations

14. Transport information

Not subject to transport regulations

15. Regulatory information

Labelling according to EC directives

Not classified as dangerous according to EC Directives

16. Other information

Issue Date: 11/10/2005

Date of print: 10/04/2008