



## TECHNICAL DATA SHEET

**DESCRIPTION** EMPTY STERILE INJECTION VIAL 10ML  
**SGL PRODUCT CODE** 0048

An empty 10ml capacity injection vial with rubber stopper and crimped flip off cap. The inside of the vial is sterilised by steam (autoclaved). There is a trace amount of sterile purified water inside the vial to facilitate steam sterilisation conditions during autoclaving.

### FORMULATION

COMPONENT
10ml Capacity Glass Injection Vial (Empty)
Rubber stopper
Flip off / Push off cap

### QUALITY CONTROL SPECIFICATION

PHYSICAL TESTS	SPECIFICATION CRITERIA
Appearance	Empty, clear and colourless glass vial with moisture visible inside

STERILITY TESTS	SPECIFICATION CRITERIA
Incubation at 20-25°C for 14 days	No growth detected
Incubation at 30-35°C for 14 days	No growth detected

Additional specification testing may be performed as requested by the customer.

Manufactured in compliance with ISO 9001 (Ref No FM37824) and tested by a UKAS (ISO 17025) accredited laboratory (Ref No. 4356).