



## TECHNICAL DATA SHEET

**DESCRIPTION** 40% GLUCOSE SOLUTION PHARMACOPOEIAL GRADE,  
BOTTLED

**SGL PRODUCT CODE** 3149

A 40% solution of glucose meeting the requirements of the European Pharmacopoeia and US Pharmacopoeia (Ph Eur/USP grade).

### FORMULATION

Typical product composition\*:

COMPONENT	WEIGHT / VOLUME
D(+) glucose, anhydrous	400 g
Purified water	600 ml

\*Product may be adjusted and/or supplemented to meet performance criteria

### QUALITY CONTROL SPECIFICATION

PHYSICAL TESTS	SPECIFICATION CRITERIA
Appearance	Clear, colourless liquid

STERILITY TESTS	SPECIFICATION CRITERIA
Incubation at 20-25°C for a minimum of 5 days	No growth detected
Incubation at 30-35°C for a minimum of 5 days	No growth detected

Additional specification testing may be performed as requested by the customer.

Manufactured in compliance with ISO 9001 (Ref No FM37824).