



TECHNICAL DATA SHEET

DESCRIPTION ISOPROPYL MYRISTATE, BOTTLED
SGL PRODUCT CODE 3377

A non-aqueous diluent used for the examination of ointments in a fatty base and emulsions of water in oil.

The product complies with the specification of the European Pharmacopoeia (EP) 2.6.13 and United States Pharmacopoeia (USP) <71>.

FORMULATION

Typical product composition*:

COMPONENT	WEIGHT / VOLUME
Isopropyl myristate	1000 ml

*Product may be adjusted and/or supplemented to meet performance criteria

QUALITY CONTROL SPECIFICATION

PHYSICAL TESTS	SPECIFICATION CRITERIA
Appearance	Clear, colourless liquid

STERILITY TESTS	SPECIFICATION CRITERIA
Incubation at 22-25°C for 5 days	No growth detected
Incubation at 35-37°C for 5 days	No growth detected
Incubation at 42-45°C for 5 days	No growth detected

Tested by addition to SGL0595 Tryptone Soya Broth

Additional specification testing may be performed as requested by the customer.

Manufactured in compliance with ISO 9001 (Ref No FM37824) and tested in accordance with ISO 11133 by a UKAS (ISO 17025) accredited laboratory (Ref No. 4356).