




CONTAINER GUIDE



CLINICAL & VETERINARY • FOOD & BEVERAGES • WATER & ENVIRONMENTAL
PHARMACEUTICAL, COSMETIC & TOILETRIES • BESPOKE MEDIA & SPECIALIST APPLICATIONS

CONVENIENCE YOU CAN TRUST

www.sglab.co.uk

Bottles, Tubes, Vials & Caps

A full range of glass bottles can be supplied with or without caps and closures in sterile or non sterile formats.

CAP KEY

Polypropylene (PP)	Standard Polypropylene Cap
Metal	Standard Metal Cap
Perforated (Perf)	Perforated Metal Cap for use with Injectable Bung

Pilfer Proof	Polypropylene Pilfer Proof Cap with Security Collar
Flip Off	Aluminium Cap with Plastic 'Flip Off' Top for use with Injectable Bung, Flip Off Top cannot be replaced.

SYRUP BOTTLES



Metal



PP



Perforated Cap, Bung and Foil



Pilfer Proof

WINCHESTER BOTTLE

With PP cap



DURAN BOTTLE



Standard Duran



Wide Neck Duran

GLASS BOTTLE GUIDE

Bottle	Material	Bottle Sizes (ml)	Neck Width	Cap Options
Syrup	Soda Glass	60, 125, 250, 300, 500, 1000	28mm	PP, Metal, Bung + Perf
Winchester	Soda Glass	500, 1000	31mm	PP
Duran Bottle	Pyrex Type 1	50, 100, 250, 500, 1000, 3500, 5000, 10000, 20000	45mm	PP or PP with Vent
Wide Neck Duran Bottle	Pyrex Type 1	250, 500, 1000, 2000, 3500, 5000, 10000, 20000	85mm	PP

UNIVERSAL VIALS



1LB HONEY JAR
With PP Cap



Thick walled Universal Vial
with PP or Metal Cap

Universal Vial with
PP or Metal Cap

BIJOU



Thick walled Bijou
with Metal Cap



Bijou with
PP Cap



Scintillation
Vial

GLASS BOTTLE GUIDE

Bottle	Material	Bottle Sizes (ml)	Neck Width	Cap Options
1LB Honey Jar	Soda Glass	225	60mm	PP
Universal Vial	Soda Glass	30	24mm	PP, Metal
Thick Walled Universal	Soda Glass	20	28mm	PP, Metal
Bijou	Soda Glass	7	18mm	PP, Metal
Round Bottom Tube	Soda Glass	15, 20	18mm	Vented Cap
Round Bottom Tube	Pyrex Type 1	10	15mm	Phenolic with Rubber Liner
Flat Bottom Tube	Soda Glass	12	18mm	PP, Metal

TUBES



Round Bottom Tube
with Phenolic Cap
and Rubber Wad



Flat Bottom
Tube with
PP Cap




Round Bottom
Tube with
Vent Cap

CAN'T SEE WHAT YOU NEED?

SGL CAN SOURCE A BOTTLE OR CONTAINER, CLOSURE AND PACKAGING TO FIT YOUR REQUIREMENT OR USE CUSTOMER SUPPLIED BOTTLES OR CONTAINERS.

PHARMACEUTICAL DIN BOTTLES

With Flip Off Cap



ENVIRONMENTAL CONCERN?

SGL OFFERS A RECYCLED OPTION WHERE BOTTLES CAN BE RETURNED AND RECYCLED AND NO CARDBOARD PACKAGING USED



GLASS BOTTLE GUIDE

Bottle	Material	Bottle Sizes (ml)	Neck Width	Cap Options
DIN Bottle	Pyrex Type 1	50, 100, 500, 1000	32mm	Flip Top with Rubber Bung
Injection Vial	Pyrex Type 1	10, 20, 50, 100	20mm	Flip Top with Rubber Bung
Plastic Stain Bottle - Square	HD	250, 1200	12mm, 20mm	PP Cap and bung
Plastic Stain Bottle - Round	HDPE	500	12mm	PP Cap and Pouring Ring
Plastic Stain Bottle - Round	Clear PVC	500	12mm	PP Cap and Pouring Ring
Plastic Chemical Bottle	HDPE	500	45mm	PP Cap and Red Safety Ring

INJECTION VIALS

With Flip Off Cap



PLASTIC BOTTLES



Round Stain Bottle

Clear Stain Bottle

Square Stain Bottle

Square Stain Bottle

Round Chemical Bottle






Plate Type	Dish Diameter	Fill Range	Print Options
55mm	55mm	→10-15ml	Plate base
55mm Contact	55mm	14ml	Plate base
90mm	90mm	→18-30ml	Bottom/side
140mm	140mm	→65-100ml	Bottom/side
Bi Plate	90mm	-	-
Tri Plate	90mm	-	-

Plate identification can be specified as side or bottom printed



Triple wrapped and irradiated plates





SGL

Southern Group Laboratory Limited

E-H Cavendish Courtyard, Sallow Road, Weldon North Industrial Estate,
Corby, Northamptonshire NN17 5JX United Kingdom

T: +44 (0)1536 403815 **F:** +44 (0)1536 403814 **E:** info@sglab.co.uk **www.sglab.co.uk**

