



TECHNICAL DATA SHEET

DESCRIPTION OXALIC ACID 3%, BOTTLED

SGL PRODUCT CODE 3109

3% oxalic acid is used for the analysis of sputum samples for *Mycobacterium* Spp..

FORMULATION

Typical product composition*:

COMPONENT	WEIGHT / VOLUME
Oxalic acid	30.0 g
Purified water	1000 ml

*Product may be adjusted and/or supplemented to meet performance criteria

QUALITY CONTROL SPECIFICATION

PHYSICAL TESTS	SPECIFICATION CRITERIA
Appearance	Clear, colourless liquid

STERILITY TESTS	SPECIFICATION CRITERIA
Incubation at 22-25°C for 5 days	No growth detected
Incubation at 35-37°C for 5 days	No growth detected
Incubation at 42-45°C for 5 days	No growth detected

Additional specification testing may be performed as requested by the customer.

Manufactured in compliance with ISO 9001 (Ref No FM37824) and tested in accordance with ISO 11133 by a UKAS (ISO 17025) accredited laboratory (Ref No. 4356).